



*tasarimmed*<sup>®</sup>

# FN-2 RECON INTRAMEDULLARY NAIL Catalogue

## DESCRIPTION

Tasarimmed FN-2 Recon Femoral nail is a cannulated intramedullary nail system with a wide range of different sized nails, Locking Screws, LAG Screws and end caps in different lengths. All implants are manufactured from ISO 5832/3 Ti6Al4V-ELI surgical alloy

## INDICATIONS

The indications of the FN-2 Recon Femoral Nail system is not different from the indications of any recon femoral nail. It is designed for use in the treatment of femur fractures in patients whose bone age is determined to be adult.

Tasarimmed Recon Femoral Nail is indicated for

- Subtrochanteric fractures
- Intertrochanteric fractures
- Ipsilateral neck / shaft fractures
- Partial proximal femur shaft fractures
- Requirement of femoral fixation after pathologic disease
- Temporary stabilization of femur shaft fractures - from the femur neck to the femur supracondylar region

## CONTRAINDICATIONS

1. These systems should not be used in crossing open epiphyseal plates.
2. Insufficient quantity or quality of bone, obliterated medullary canal or conditions which tend to retard healing, blood supply limitations, previous infections, etc.
3. Active infection.
4. Any hardware that would preclude use of nails.
5. Congenital or acquired bony deformity.
6. Hypovolemia, hypothermia and coagulopathy.
7. Mental conditions that preclude cooperation with the rehabilitation regimen.
8. Morbid obesity
9. Alcoholism, smoking or drug abuse and addiction
10. Fever or leukocytosis.
11. Pregnancy.
12. Metal sensitivities / allergies.

13. Severe osteoporosis.
14. Congenital abnormalities, tumors or the presence of any condition that prevents safe implantation and / or shortens the life span.
15. Any condition, including inadequate tissue coverage on the operation site where the device will interfere with the anatomical structure or physiological performance.
16. Patients who do not want or cannot follow post-operative care instructions.
17. Open fractures
18. Any situation that is not identified under the heading INDICATIONS.

## Implant List


Info	Ref No	Description
Single Use	FN2L-100280	FN-2 Recon Nail LEFT Ø 10 L280 mm
Single Use	FN2L-100300	FN-2 Recon Nail LEFT Ø 10 L300 mm
Single Use	FN2L-100320	FN-2 Recon Nail LEFT Ø 10 L320 mm
Single Use	FN2L-100340	FN-2 Recon Nail LEFT Ø 10 L340 mm
Single Use	FN2L-100360	FN-2 Recon Nail LEFT Ø 10 L360 mm
Single Use	FN2L-100380	FN-2 Recon Nail LEFT Ø 10 L380 mm
Single Use	FN2L-100400	FN-2 Recon Nail LEFT Ø 10 L400 mm
Single Use	FN2L-100420	FN-2 Recon Nail LEFT Ø 10 L420 mm
Single Use	FN2L-110280	FN-2 Recon Nail LEFT Ø 11 L280 mm
Single Use	FN2L-110300	FN-2 Recon Nail LEFT Ø 11 L300 mm
Single Use	FN2L-110320	FN-2 Recon Nail LEFT Ø 11 L320 mm
Single Use	FN2L-110340	FN-2 Recon Nail LEFT Ø 11 L340 mm
Single Use	FN2L-110360	FN-2 Recon Nail LEFT Ø 11 L360 mm
Single Use	FN2L-110380	FN-2 Recon Nail LEFT Ø 11 L380 mm
Single Use	FN2L-110400	FN-2 Recon Nail LEFT Ø 11 L400 mm
Single Use	FN2L-110420	FN-2 Recon Nail LEFT Ø 11 L420 mm
Single Use	FN2L-125280	FN-2 Recon Nail LEFT Ø 12.5 L280 mm
Single Use	FN2L-125300	FN-2 Recon Nail LEFT Ø 12.5 L300 mm
Single Use	FN2L-125320	FN-2 Recon Nail LEFT Ø 12.5 L320 mm
Single Use	FN2L-125340	FN-2 Recon Nail LEFT Ø 12.5 L340 mm
Single Use	FN2L-125360	FN-2 Recon Nail LEFT Ø 12.5 L360 mm
Single Use	FN2L-125380	FN-2 Recon Nail LEFT Ø 12.5 L380 mm
Single Use	FN2L-125400	FN-2 Recon Nail LEFT Ø 12.5 L400 mm
Single Use	FN2L-125420	FN-2 Recon Nail LEFT Ø 12.5 L420 mm
Single Use	FN2L-140280	FN-2 Recon Nail LEFT Ø 14 L280 mm
Single Use	FN2L-140300	FN-2 Recon Nail LEFT Ø 14 L300 mm
Single Use	FN2L-140320	FN-2 Recon Nail LEFT Ø 14 L320 mm
Single Use	FN2L-140340	FN-2 Recon Nail LEFT Ø 14 L340 mm
Single Use	FN2L-140360	FN-2 Recon Nail LEFT Ø 14 L360 mm
Single Use	FN2L-140380	FN-2 Recon Nail LEFT Ø 14 L380 mm
Single Use	FN2L-140400	FN-2 Recon Nail LEFT Ø 14 L400 mm
Single Use	FN2L-140420	FN-2 Recon Nail LEFT Ø 14 L420 mm



Info	Ref No	Description
Single Use	FN2R-100280	FN-2 Recon Nail RIGHT Ø 10 L280 mm
Single Use	FN2R-100300	FN-2 Recon Nail RIGHT Ø 10 L300 mm
Single Use	FN2R-100320	FN-2 Recon Nail RIGHT Ø 10 L320 mm
Single Use	FN2R-100340	FN-2 Recon Nail RIGHT Ø 10 L340 mm
Single Use	FN2R-100360	FN-2 Recon Nail RIGHT Ø 10 L360 mm
Single Use	FN2R-100380	FN-2 Recon Nail RIGHT Ø 10 L380 mm
Single Use	FN2R-100400	FN-2 Recon Nail RIGHT Ø 10 L400 mm
Single Use	FN2R-100420	FN-2 Recon Nail RIGHT Ø 10 L420 mm
Single Use	FN2R-110280	FN-2 Recon Nail RIGHT Ø 11 L280 mm
Single Use	FN2R-110300	FN-2 Recon Nail RIGHT Ø 11 L300 mm
Single Use	FN2R-110320	FN-2 Recon Nail RIGHT Ø 11 L320 mm
Single Use	FN2R-110340	FN-2 Recon Nail RIGHT Ø 11 L340 mm
Single Use	FN2R-110360	FN-2 Recon Nail RIGHT Ø 11 L360 mm
Single Use	FN2R-110380	FN-2 Recon Nail RIGHT Ø 11 L380 mm
Single Use	FN2R-110400	FN-2 Recon Nail RIGHT Ø 11 L400 mm
Single Use	FN2R-110420	FN-2 Recon Nail RIGHT Ø 11 L420 mm
Single Use	FN2R-125280	FN-2 Recon Nail RIGHT Ø 12.5 L280 mm
Single Use	FN2R-125300	FN-2 Recon Nail RIGHT Ø 12.5 L300 mm
Single Use	FN2R-125320	FN-2 Recon Nail RIGHT Ø 12.5 L320 mm
Single Use	FN2R-125340	FN-2 Recon Nail RIGHT Ø 12.5 L340 mm
Single Use	FN2R-125360	FN-2 Recon Nail RIGHT Ø 12.5 L360 mm
Single Use	FN2R-125380	FN-2 Recon Nail RIGHT Ø 12.5 L380 mm
Single Use	FN2R-125400	FN-2 Recon Nail RIGHT Ø 12.5 L400 mm
Single Use	FN2R-125420	FN-2 Recon Nail RIGHT Ø 12.5 L420 mm
Single Use	FN2R-140280	FN-2 Recon Nail RIGHT Ø 14 L280 mm
Single Use	FN2R-140300	FN-2 Recon Nail RIGHT Ø 14 L300 mm
Single Use	FN2R-140320	FN-2 Recon Nail RIGHT Ø 14 L320 mm
Single Use	FN2R-140340	FN-2 Recon Nail RIGHT Ø 14 L340 mm
Single Use	FN2R-140360	FN-2 Recon Nail RIGHT Ø 14 L360 mm
Single Use	FN2R-140380	FN-2 Recon Nail RIGHT Ø 14 L380 mm
Single Use	FN2R-140400	FN-2 Recon Nail RIGHT Ø 14 L400 mm
Single Use	FN2R-140420	FN-2 Recon Nail RIGHT Ø 14 L420 mm













Info	Ref No	Description
Single Use	FN2LS-7060	FN-2 LAG Screw Ø 7 L60 mm
Single Use	FN2LS-7065	FN-2 LAG Screw Ø 7 L65 mm
Single Use	FN2LS-7070	FN-2 LAG Screw Ø 7 L70 mm
Single Use	FN2LS-7075	FN-2 LAG Screw Ø 7 L75 mm
Single Use	FN2LS-7080	FN-2 LAG Screw Ø 7 L80 mm
Single Use	FN2LS-7085	FN-2 LAG Screw Ø 7 L85 mm
Single Use	FN2LS-7090	FN-2 LAG Screw Ø 7 L90 mm
Single Use	FN2LS-7095	FN-2 LAG Screw Ø 7 L95 mm
Single Use	FN2LS-7100	FN-2 LAG Screw Ø 7 L100 mm
Single Use	FN2LS-7105	FN-2 LAG Screw Ø 7 L105 mm
Single Use	FN2LS-7110	FN-2 LAG Screw Ø 7 L110 mm
Single Use	FN2LS-7115	FN-2 LAG Screw Ø 7 L115 mm
Single Use	FN2LS-7120	FN-2 LAG Screw Ø 7 L120 mm

End Cap		
		
Info	Ref No	Description
Single Use	FIN-1001	FN-2 Nail End Cap 1 mm
Single Use	FIN-1005	FN-2 Nail End Cap 5 mm
Single Use	FIN-1010	FN-2 Nail End Cap 10 mm
Single Use	FIN-1015	FN-2 Nail End Cap 15 mm
Single Use	FIN-1020	FN-2 Nail End Cap 20 mm



Info	Ref No	Description
Single Use	INS-5024F	Nail Locking Screw Ø5 mm L24 (Full Threaded)
Single Use	INS-5026F	Nail Locking Screw Ø5 mm L26 (Full Threaded)
Single Use	INS-5028F	Nail Locking Screw Ø5 mm L28 (Full Threaded)
Single Use	INS-5030F	Nail Locking Screw Ø5 mm L30 (Full Threaded)
Single Use	INS-5032F	Nail Locking Screw Ø5 mm L32 (Full Threaded)
Single Use	INS-5034F	Nail Locking Screw Ø5 mm L34 (Full Threaded)
Single Use	INS-5036F	Nail Locking Screw Ø5 mm L36 (Full Threaded)
Single Use	INS-5038F	Nail Locking Screw Ø5 mm L38 (Full Threaded)
Single Use	INS-5040F	Nail Locking Screw Ø5 mm L40 (Full Threaded)
Single Use	INS-5045F	Nail Locking Screw Ø5 mm L45 (Full Threaded)
Single Use	INS-5050F	Nail Locking Screw Ø5 mm L50 (Full Threaded)
Single Use	INS-5055F	Nail Locking Screw Ø5 mm L55 (Full Threaded)
Single Use	INS-5060F	Nail Locking Screw Ø5 mm L60 (Full Threaded)
Single Use	INS-5065F	Nail Locking Screw Ø5 mm L65 (Full Threaded)
Single Use	INS-5070F	Nail Locking Screw Ø5 mm L70 (Full Threaded)
Single Use	INS-5075F	Nail Locking Screw Ø5 mm L75 (Full Threaded)
Single Use	INS-5080F	Nail Locking Screw Ø5 mm L80 (Full Threaded)
Single Use	INS-5085F	Nail Locking Screw Ø5 mm L85 (Full Threaded)
Single Use	INS-5085F	Nail Locking Screw Ø5 mm L85 (Full Threaded)
Single Use	INS-5090F	Nail Locking Screw Ø5 mm L90 (Full Threaded)
Single Use	INS-5095F	Nail Locking Screw Ø5 mm L95 (Full Threaded)
Single Use	INS-5100F	Nail Locking Screw Ø5 mm L100 (Full Threaded)
Single Use	INS-5105F	Nail Locking Screw Ø5 mm L105 (Full Threaded)
Single Use	INS-5110F	Nail Locking Screw Ø5 mm L110 (Full Threaded)

## Symbols used in Packaging and their meanings

-  Indicates a medical device that has not been subjected to a sterilization process.
-  Indicates the manufacturer's batch code so that the batch or lot can be identified  
NOTE: Synonyms for "batch code" are "lot number" and "batch number"
-  Indicates the manufacturer's catalog number so that the medical device can be identified  
NOTE: Synonyms for "catalogue number" are "reference number" and "reorder number"
-  Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
-  Indicates the need for the user to consult the instructions for use.
-  Indicates a medical device that should not be used if the package has been damaged or opened.
-  Indicates the medical device manufacturer, as defined in EU Directives
-  Symbol for date of manufacture. This symbol is accompanied by a date.
-  Indicates the date after which the medical device is not to be used.
-  Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

---

B04-TDTG03 FN-2 Nail Cat EN

All rights reserved. Tasarimmed FN-2 Product Catalogue Rev.0 16.07.2020

Tasarimmed Tibbi Mamuller San. Tic. A.S.

CE1984

Topcular C. Setustu no.3 34055 Eyup, Istanbul Turkey [www.tasarimmed.com.tr](http://www.tasarimmed.com.tr)

Tel: 0090 212 674 22 44-45 Fax: 0090 212 674 2246

---